

## Elranatamab

## Chemical Properties

CAS No. : 2408850-14-4

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Elranatamab (PF-06863135) is a humanized bispecific antibody that targets BCMA (B-cell maturation antigen) on myeloma cells and CD3 on T cells. By targeting these molecules, Elranatamab activates T cells and directs them to induce a cytotoxic T cell response against myeloma cells. This compound is used in research for relapsed or refractory multiple myeloma.
Targets(IC50)	TNF

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481