

HH0043

Chemical Properties


CAS No. : 2922069-42-7

Formula: C₂₃H₂₄F₃N₅O

Molecular Weight: 443.46

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	HH0043 (Compound 10f) is an orally active and potent inhibitor of SOS1, exhibiting an IC ₅₀ of 5.8 nM. It is applicable in cancer research.
Targets(IC ₅₀)	Ras

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.255 mL	11.275 mL	22.5499 mL
5 mM	0.451 mL	2.255 mL	4.510 mL
10 mM	0.2255 mL	1.1275 mL	2.255 mL
50 mM	0.0451 mL	0.2255 mL	0.451 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481