


Senkyunolide E

Chemical Properties

| | |
|-------------------|--|
| CAS No. : | 94530-83-3 |
| Formula: | C ₁₂ H ₁₂ O ₃ |
| Molecular Weight: | 204.22 |
| Storage: | Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small> |



Biological Description

| | |
|-------------|---|
| Description | Senkyunolide E is a natural product for research related to life sciences. The catalog number is TN5818 and the CAS number is 94530-83-3. |
|-------------|---|

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 4.8967 mL | 24.4834 mL | 48.9668 mL |
| 5 mM | 0.9793 mL | 4.8967 mL | 9.7934 mL |
| 10 mM | 0.4897 mL | 2.4483 mL | 4.8967 mL |
| 50 mM | 0.0979 mL | 0.4897 mL | 0.9793 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481