

Azido-PEG5-CH<sub>2</sub>CO<sub>2</sub>-NHS

## Chemical Properties

CAS No. : 2144777-77-3


Formula: C<sub>16</sub>H<sub>26</sub>N<sub>4</sub>O<sub>9</sub>

Molecular Weight: 418.4

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Azido-PEG5-CH <sub>2</sub> CO <sub>2</sub> -NHS is a PEG-based linker for PROTACs that bonds two essential ligands, facilitating the formation of PROTAC molecules and enabling selective protein degradation through the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs, which consist of two ligands connected by a linker—one binding to an E3 ubiquitin ligase and the other to the target protein—utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3901 mL	11.9503 mL	23.9006 mL
5 mM	0.478 mL	2.3901 mL	4.7801 mL
10 mM	0.239 mL	1.195 mL	2.3901 mL
50 mM	0.0478 mL	0.239 mL	0.478 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481