

(2S,3S,4S,5R,6S)-3,4,5-trihydroxy-6-({4-phenyl-[1,1'-biphenyl]-2,6-diyl}oxy)oxane-2-carboxylic acid


## Chemical Properties

CAS No. :

Formula: C<sub>24</sub>H<sub>22</sub>O<sub>9</sub>

Molecular Weight: 454.43

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(2S,3S,4S,5R,6S)-3,4,5-trihydroxy-6-({4-phenyl-[1,1'-biphenyl]-2,6-diyl}oxy)oxane-2-carboxylic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2006 mL	11.0028 mL	22.0056 mL
5 mM	0.4401 mL	2.2006 mL	4.4011 mL
10 mM	0.2201 mL	1.1003 mL	2.2006 mL
50 mM	0.044 mL	0.2201 mL	0.4401 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481