

Spectabiline

Chemical Properties


CAS No. :

Formula: C₁₈H₂₅N₀O₇

Molecular Weight: 367.398

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Spectabiline is a useful organic compound for research related to life sciences and the catalog number is T124728.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7218 mL	13.6091 mL	27.2183 mL
5 mM	0.5444 mL	2.7218 mL	5.4437 mL
10 mM	0.2722 mL	1.3609 mL	2.7218 mL
50 mM	0.0544 mL	0.2722 mL	0.5444 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481