


1,2,3-Triphenylguanidine

Chemical Properties

CAS No. :	101-01-9
Formula:	C ₁₉ H ₁₇ N ₃
Molecular Weight:	287.37
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	1,2,3-Triphenylguanidine is an agent of biochemical.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4798 mL	17.3992 mL	34.7983 mL
5 mM	0.696 mL	3.4798 mL	6.9597 mL
10 mM	0.348 mL	1.7399 mL	3.4798 mL
50 mM	0.0696 mL	0.348 mL	0.696 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481