


## 1-Benzyl-4-hydroxypiperidine

## Chemical Properties

CAS No. :	4727-72-4
Formula:	C <sub>12</sub> H <sub>17</sub> NO
Molecular Weight:	191.274
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1-Benzyl-4-hydroxypiperidine is used as a pharmaceutical intermediate for the synthesis of many drugs.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2282 mL	26.1411 mL	52.2821 mL
5 mM	1.0456 mL	5.2282 mL	10.4564 mL
10 mM	0.5228 mL	2.6141 mL	5.2282 mL
50 mM	0.1046 mL	0.5228 mL	1.0456 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481