

GST Tag Agarose

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at 4°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Targetmol's GST Tag Agarose use a highly cross-linked agarose gel as a matrix, chemically modified with a spacer arm of 12 atoms in length to covalently bind reduced glutathione. This design ensures high protein binding capacity and stability, and GST Tag Agarose Gels are specifically designed for the purification of GST Tag fusion recombinant proteins from different expression systems.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481