

Flag Tag Magrose Beads

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at 4°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Targetmol's Flag Tag Magrose Beads covalently bind a mouse-derived anti-DYKDDDDK (Flag) monoclonal antibody to the surface of the agarose beads, which specifically binds to Flag-tagged proteins.

Targetmol's Flag Tag Magrose Beads can be used for immunoprecipitation (IP), co-immunoprecipitation (Co-IP) or Pull-down assays of proteins and protein complexes, as well as for purification of Flag tagged fusion proteins. It is suitable for antigenic samples from cell lysates, cell secretion supernatants, serum, ascites and other sources.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481