

GST Tag Magrose Beads

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at 4°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Targetmol's GST Tag Magrose Beads use highly cross-linked agarose magnetic beads as a matrix, chemically modified with a 12-atom length spacer arm to covalently bind reduced glutathione. This design ensures efficient protein binding and stability. GST Tag Agarose Beads can be used to purify GST-tagged fusion recombinant proteins from different expression systems, as well as for Pull-down assays of GST-tagged proteins.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481