


## 20S,24R-Epoxydammar-12,25-diol-3-one

### Chemical Properties

CAS No. :	25279-15-6
Formula:	C <sub>30</sub> H <sub>50</sub> O <sub>4</sub>
Molecular Weight:	474.726
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



### Biological Description

Description	20S,24R-Epoxydammar-12,25-diol-3-one is a natural product of Pseudolarix, Pinaceae. The catalog number is TN2803 and the CAS number is 25279-15-6. 20S,24R-Epoxydammar-12,25-diol-3-one can be used as a reference standard.
Targets(IC50)	Others

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1065 mL	10.5323 mL	21.0646 mL
5 mM	0.4213 mL	2.1065 mL	4.2129 mL
10 mM	0.2106 mL	1.0532 mL	2.1065 mL
50 mM	0.0421 mL	0.2106 mL	0.4213 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481