


## SMARCA2 ligand-6

## Chemical Properties

CAS No. :	2408394-89-6
Formula:	C <sub>18</sub> H <sub>20</sub> N <sub>6</sub> O
Molecular Weight:	336.39
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	SMARCA2 ligand-6 is a PROTAC-targeted protein ligand.
-------------	-------------------------------------------------------

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9727 mL	14.8637 mL	29.7274 mL
5 mM	0.5945 mL	2.9727 mL	5.9455 mL
10 mM	0.2973 mL	1.4864 mL	2.9727 mL
50 mM	0.0595 mL	0.2973 mL	0.5945 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481