

## COOEt-spiro[3.3]heptane-CHO

## Chemical Properties

CAS No. : 2226571-47-5

Formula: C<sub>11</sub>H<sub>16</sub>O<sub>3</sub>

Molecular Weight: 196.24

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	PROTACSMARCA2 degrader-6-Linker is a PROTAC linker that joins an E3 ligase ligand with a SMARCA2 ligand to form the PROTACSMARCA2 degrader-6. This compound enables the targeted degradation of the SMARCA2 protein by recruiting an E3 ligase for ubiquitination and subsequent proteasomal degradation.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.0958 mL	25.479 mL	50.958 mL
5 mM	1.0192 mL	5.0958 mL	10.1916 mL
10 mM	0.5096 mL	2.5479 mL	5.0958 mL
50 mM	0.1019 mL	0.5096 mL	1.0192 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481