


Abenacianine

Chemical Properties

CAS No. :	2231255-31-3
Formula:	C127H145ClF4N10O23S5
Molecular Weight:	2451.33
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Abenacianine is a diagnostic imaging agent.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.4079 mL	2.0397 mL	4.0794 mL
5 mM	0.0816 mL	0.4079 mL	0.8159 mL
10 mM	0.0408 mL	0.204 mL	0.4079 mL
50 mM	0.0082 mL	0.0408 mL	0.0816 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481