

Cemavafusp

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cemavafusp is a monoclonal antibody targeting PD-L1, composed of a single chain variable fragment (single chain variable fragment), an antibody heavy chain, and a κ light chain.
Targets(IC50)	PD-1/PD-L1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481