

SMARCA2/4-ligand-2

Chemical Properties

CAS No. : 2568281-88-7


Formula: C₂₀H₂₀BrN₇O

Molecular Weight: 454.32

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

SMARCA2/4-ligand-2 is a PROTAC target protein ligand.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2011 mL	11.0055 mL	22.0109 mL
5 mM	0.4402 mL	2.2011 mL	4.4022 mL
10 mM	0.2201 mL	1.1005 mL	2.2011 mL
50 mM	0.044 mL	0.2201 mL	0.4402 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481