

HKYK-0030

## Chemical Properties


CAS No. : 2923054-94-6

Formula: C<sub>14</sub>H<sub>11</sub>ClN<sub>4</sub>O<sub>5</sub>S

Molecular Weight: 382.78

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	HKYK-0030 is an inhibitor of the MTDH-SND1 protein interaction and possesses activity that inhibits cancer metastasis.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6125 mL	13.0623 mL	26.1247 mL
5 mM	0.5225 mL	2.6125 mL	5.2249 mL
10 mM	0.2612 mL	1.3062 mL	2.6125 mL
50 mM	0.0522 mL	0.2612 mL	0.5225 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481