

Alloc-Val-Ala-PAB-PNP

Chemical Properties


CAS No. : 1884578-27-1

Formula: C₂₆H₃₀N₄O₉

Molecular Weight: 542.54

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Alloc-Val-Ala-PAB-PNP is a cleavable ADC linker.
Targets(IC50)	ADC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8432 mL	9.2159 mL	18.4318 mL
5 mM	0.3686 mL	1.8432 mL	3.6864 mL
10 mM	0.1843 mL	0.9216 mL	1.8432 mL
50 mM	0.0369 mL	0.1843 mL	0.3686 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481