


Segetalin C

Chemical Properties

CAS No. :	164991-90-6
Formula:	C40H51N9O7
Molecular Weight:	769.89
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Segetalin C, a cyclic peptide, effectively increases uterine weight in ovariectomized rats (2.5 mg/kg).
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2989 mL	6.4944 mL	12.9889 mL
5 mM	0.2598 mL	1.2989 mL	2.5978 mL
10 mM	0.1299 mL	0.6494 mL	1.2989 mL
50 mM	0.026 mL	0.1299 mL	0.2598 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481