

JP-2-249

## Chemical Properties

CAS No. : 3036530-41-0


Formula: C<sub>25</sub>H<sub>25</sub>N<sub>5</sub>O<sub>4</sub>

Molecular Weight: 459.5

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	JP-2-249, a molecular glue, serves as an effective degrader of SMARCA2. At concentrations ranging from 1-10 $\mu$ M, JP-2-249 reduces the protein levels of SMARCA2 in MV-4-11 cells.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1763 mL	10.8814 mL	21.7628 mL
5 mM	0.4353 mL	2.1763 mL	4.3526 mL
10 mM	0.2176 mL	1.0881 mL	2.1763 mL
50 mM	0.0435 mL	0.2176 mL	0.4353 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481