


PROTAC SMARCA2/4-degrader-11

Chemical Properties

CAS No. :	2568276-71-9
Formula:	C37H41N9O4
Molecular Weight:	675.78
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	SMARCA2 degrader-1 (compound I-321) is a PROTAC degrader with a DC50 for SMARCA of less than 100nM, featuring a structure that incorporates a SMARCA inhibitor, a linker, and an E3 ligase ligand.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4798 mL	7.3989 mL	14.7977 mL
5 mM	0.296 mL	1.4798 mL	2.9595 mL
10 mM	0.148 mL	0.7399 mL	1.4798 mL
50 mM	0.0296 mL	0.148 mL	0.296 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481