

Dolichol-20

Chemical Properties

CAS No. :	2067-66-5
Formula:	C100H164O
Molecular Weight:	1382.37
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Dolichol-20 is categorized as an ester product.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7234 mL	3.617 mL	7.234 mL
5 mM	0.1447 mL	0.7234 mL	1.4468 mL
10 mM	0.0723 mL	0.3617 mL	0.7234 mL
50 mM	0.0145 mL	0.0723 mL	0.1447 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481