


## PROTAC SMARCA2 degrader-8

## Chemical Properties

CAS No. :	2669027-50-1
Formula:	C43H50N10O5S
Molecular Weight:	818.99
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	PROTACSMARCA2 degrader-8 (compound 1) is an effective SMARCA2 degrader, exhibiting a DC50 of 28nM in A375 cells.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.221 mL	6.1051 mL	12.2102 mL
5 mM	0.2442 mL	1.221 mL	2.442 mL
10 mM	0.1221 mL	0.6105 mL	1.221 mL
50 mM	0.0244 mL	0.1221 mL	0.2442 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481