


SMARCA2/4-ligand-1

Chemical Properties

CAS No. :	2411651-36-8
Formula:	C ₁₉ H ₂₂ N ₆ O
Molecular Weight:	350.42
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	SMARCA2/4-ligand-1 is a PROTAC target protein ligand.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8537 mL	14.2686 mL	28.5372 mL
5 mM	0.5707 mL	2.8537 mL	5.7074 mL
10 mM	0.2854 mL	1.4269 mL	2.8537 mL
50 mM	0.0571 mL	0.2854 mL	0.5707 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481