

SMARCA2 degrader-20

Chemical Properties

CAS No. : 2568273-85-6


Formula: C₃₈H₄₁N₉O₅

Molecular Weight: 703.79

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

SMARCA2 degrader-20, a potent PROTAC SMARCA2 degrader, exhibits a DC50 value of less than 100 nM in A549 cells.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4209 mL	7.1044 mL	14.2088 mL
5 mM	0.2842 mL	1.4209 mL	2.8418 mL
10 mM	0.1421 mL	0.7104 mL	1.4209 mL
50 mM	0.0284 mL	0.1421 mL	0.2842 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481