

SMARCA2 ligand-5

Chemical Properties

CAS No. : 2408394-16-9


Formula: C₂₅H₃₀N₈O

Molecular Weight: 458.56

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

SMARCA2 ligand-5 is a PROTAC (proteolysis targeting chimera) target protein ligand.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1807 mL	10.9037 mL	21.8074 mL
5 mM	0.4361 mL	2.1807 mL	4.3615 mL
10 mM	0.2181 mL	1.0904 mL	2.1807 mL
50 mM	0.0436 mL	0.2181 mL	0.4361 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481