


L-Hydroxy arginine dihydrochloride

Chemical Properties

CAS No. :	133565-48-7
Formula:	C ₆ H ₁₆ Cl ₂ N ₄ O ₃
Molecular Weight:	263.12
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <i>Actual storage temperature shall be subject to the COA.</i>



Biological Description

Description	L-Hydroxy arginine dihydrochloride is a stable intermediate in the decomposition metabolism of L-arginine.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8005 mL	19.0027 mL	38.0055 mL
5 mM	0.7601 mL	3.8005 mL	7.6011 mL
10 mM	0.3801 mL	1.9003 mL	3.8005 mL
50 mM	0.076 mL	0.3801 mL	0.7601 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481