


Talogreptide mesaroxetan

Chemical Properties

CAS No. :	1801418-23-4
Formula:	C ₈₆ H ₁₄₀ N ₂₂ O ₁₈
Molecular Weight:	1770.17
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Talogreptide mesaroxetan is a diagnostic imaging agent that exhibits antitumor activity.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.5649 mL	2.8246 mL	5.6492 mL
5 mM	0.113 mL	0.5649 mL	1.1298 mL
10 mM	0.0565 mL	0.2825 mL	0.5649 mL
50 mM	0.0113 mL	0.0565 mL	0.113 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481