

Pomalidomide 5-piperidylamine

Chemical Properties

CAS No. : 2636798-55-3


Formula: C₁₈H₂₀N₄O₄

Molecular Weight: 356.38

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pomalidomide 5-piperidylamine serves as a conjugate of a linker and an E3 ligase ligand. This compound is utilized in the synthesis of PROTACPARP1 degrader-2.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.806 mL	14.030 mL	28.0599 mL
5 mM	0.5612 mL	2.806 mL	5.612 mL
10 mM	0.2806 mL	1.403 mL	2.806 mL
50 mM	0.0561 mL	0.2806 mL	0.5612 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481