


## Icaritin-3-O-rhamnoside

## Chemical Properties

CAS No. :	108195-76-2
Formula:	C <sub>27</sub> H <sub>32</sub> O <sub>11</sub>
Molecular Weight:	532.54
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Icaritin-3-O-rhamnoside (compound 20) is a flavonoid compound that exhibits DPPH free radical scavenging activity (IC <sub>50</sub> =227.9 μM).
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8778 mL	9.389 mL	18.7779 mL
5 mM	0.3756 mL	1.8778 mL	3.7556 mL
10 mM	0.1878 mL	0.9389 mL	1.8778 mL
50 mM	0.0376 mL	0.1878 mL	0.3756 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481