


## Z-TPE-2DEVD

## Chemical Properties

CAS No. :	1426396-40-8
Formula:	C78H94N18O24
Molecular Weight:	1667.69
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Z-TPE-2DEVD is a specific probe for caspase-3.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.5996 mL	2.9982 mL	5.9963 mL
5 mM	0.1199 mL	0.5996 mL	1.1993 mL
10 mM	0.060 mL	0.2998 mL	0.5996 mL
50 mM	0.012 mL	0.060 mL	0.1199 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481