


Arg-Phe-Asp-Ser

Chemical Properties

CAS No. :	102567-19-1
Formula:	C ₂₂ H ₃₃ N ₇ O ₈
Molecular Weight:	523.54
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Arg-Phe-Asp-Ser is a cleavable ADC linker.
Targets(IC50)	ADC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9101 mL	9.5504 mL	19.1007 mL
5 mM	0.382 mL	1.9101 mL	3.8201 mL
10 mM	0.191 mL	0.955 mL	1.9101 mL
50 mM	0.0382 mL	0.191 mL	0.382 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481