


Exatecan Intermediate 5

Chemical Properties

CAS No. :	143655-70-3
Formula:	C ₁₅ H ₁₇ FN ₂ O ₃
Molecular Weight:	292.31
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Exatecan Intermediate 5 serves as an intermediate in the synthesis of Exatecan (DX-8951), a DNA topoisomerase I (topoisomerase I) inhibitor with an IC ₅₀ of 2.2 μM (0.975 μg/mL), utilized in cancer research. This intermediate is also employed in the production of Antibody-Drug Conjugates (ADCs).
Targets(IC ₅₀)	ADC Cytotoxin

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.421 mL	17.1051 mL	34.2103 mL
5 mM	0.6842 mL	3.421 mL	6.8421 mL
10 mM	0.3421 mL	1.7105 mL	3.421 mL
50 mM	0.0684 mL	0.3421 mL	0.6842 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481