


## N3-L-Lys(Mtt)-OH

## Chemical Properties

CAS No. :	1333231-26-7
Formula:	C <sub>26</sub> H <sub>28</sub> N <sub>4</sub> O <sub>2</sub>
Molecular Weight:	428.53
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	N3-L-Lys(Mtt)-OH is a click chemistry reagent featuring an azido group. It efficiently introduces an azido functional group into the N-terminal position of amino acid structural units within peptide sequences.
Targets(IC50)	ADC Linker

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3336 mL	11.6678 mL	23.3356 mL
5 mM	0.4667 mL	2.3336 mL	4.6671 mL
10 mM	0.2334 mL	1.1668 mL	2.3336 mL
50 mM	0.0467 mL	0.2334 mL	0.4667 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481