

Fmoc-L-Asn(EDA-N3)-OH

Chemical Properties


CAS No. : 2616562-74-2

Formula: C₂₁H₂₁N₅O₅

Molecular Weight: 423.42

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Fmoc-L-Asn(EDA-N3)-OH is a chemical compound equipped with an azide group used in click chemistry. Reports in the literature have demonstrated its application in the modification of α -amanitin, the deadliest member of the amatoxin peptide family produced by the death cap mushroom. α -Amanitin is a highly lethal natural compound (LD ₅₀ =50-100 μ g/kg), an orally ingested rigid bicyclic octapeptide, acting as a highly selective conformational inhibitor of RNA polymerase II.
Targets(IC50)	ADC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3617 mL	11.8086 mL	23.6172 mL
5 mM	0.4723 mL	2.3617 mL	4.7234 mL
10 mM	0.2362 mL	1.1809 mL	2.3617 mL
50 mM	0.0472 mL	0.2362 mL	0.4723 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481