

## Boc-L-Lys(N3)-OH CHA

## Chemical Properties


CAS No. : 2098497-30-2

Formula: C17H33N5O4

Molecular Weight: 371.48

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Boc-L-Lys(N3)-OH CHA is a chemical reagent characterized by the presence of an azide group.
Targets(IC50)	ADC Linker

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6919 mL	13.4597 mL	26.9193 mL
5 mM	0.5384 mL	2.6919 mL	5.3839 mL
10 mM	0.2692 mL	1.346 mL	2.6919 mL
50 mM	0.0538 mL	0.2692 mL	0.5384 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481