


## Azido-PEG(4)-Val-Cit-PAB-PNP

## Chemical Properties

CAS No. :	1869126-60-2
Formula:	C <sub>36</sub> H <sub>51</sub> N <sub>9</sub> O <sub>13</sub>
Molecular Weight:	817.84
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Azido-PEG(4)-Val-Cit-PAB-PNP is a click chemistry reagent that incorporates an azide group.
Targets(IC50)	ADC Linker

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2227 mL	6.1137 mL	12.2273 mL
5 mM	0.2445 mL	1.2227 mL	2.4455 mL
10 mM	0.1223 mL	0.6114 mL	1.2227 mL
50 mM	0.0245 mL	0.1223 mL	0.2445 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481