

N3-Gly-Gly-Gly-Gly-Gly-OH

Chemical Properties


CAS No. : 2250433-77-1

Formula: C10H15N7O6

Molecular Weight: 329.27

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	N3-Gly-Gly-Gly-Gly-Gly-OH is a click chemistry reagent featuring an azido group.
Targets(IC50)	ADC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.037 mL	15.1851 mL	30.3702 mL
5 mM	0.6074 mL	3.037 mL	6.074 mL
10 mM	0.3037 mL	1.5185 mL	3.037 mL
50 mM	0.0607 mL	0.3037 mL	0.6074 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481