


Alkyne-SNAP

Chemical Properties

CAS No. :	1104822-07-2
Formula:	C ₁₈ H ₁₈ N ₆ O ₂
Molecular Weight:	350.38
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Alkyne-SNAP (compound 3) is an alkyne-coupled benzylguanine that interacts with SNAP-tag for irreversible covalent labeling of SNAP fusion proteins. The alkyne segment of the molecule can subsequently undergo click chemistry for further labeling.
Targets(IC50)	ADC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.854 mL	14.2702 mL	28.5404 mL
5 mM	0.5708 mL	2.854 mL	5.7081 mL
10 mM	0.2854 mL	1.427 mL	2.854 mL
50 mM	0.0571 mL	0.2854 mL	0.5708 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481