

## 10-Undecynoyl-OSu

## Chemical Properties


CAS No. : 1006592-57-9

Formula: C<sub>15</sub>H<sub>21</sub>NO<sub>4</sub>

Molecular Weight: 279.33

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	10-Undecynoyl-OSu is a click chemistry reagent utilized as a hydrophobic bioconjugation agent, allowing for further modification on the alkyne using Click-chemistry.
Targets(IC50)	ADC Linker

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.580 mL	17.900 mL	35.7999 mL
5 mM	0.716 mL	3.580 mL	7.160 mL
10 mM	0.358 mL	1.790 mL	3.580 mL
50 mM	0.0716 mL	0.358 mL	0.716 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481