

FLAG peptide TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	FLAG peptide TFA is the trifluoroacetate (TFA) salt form of FLAG peptide. It serves as a versatile fusion tag for the purification of recombinant proteins. FLAG peptide TFA allows its fusion proteins to maintain their natural folding. It can be removed by enterokinase and eluted under non-denaturing conditions.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481