

HBpep-SA acetate

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	HBpep-SA acetate is a cell membrane-permeable peptide aggregate that undergoes phase separation and forms stable droplets in acidic environments with pH below 6.5, enabling efficient delivery of biomolecules (such as therapeutic peptides and enzymes) and serving as a carrier for mRNA-based vaccines.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481