

DSPE-PEG2000-Cy5.5

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DSPE-PEG2000-Cy5.5 is a Cy5.5-labeled DSPE-PEG2000 conjugate suitable for drug delivery and fluorescence imaging. Cy5.5 is a near-infrared (NIR) fluorescent dye commonly used for in vivo imaging.
Targets(IC50)	Liposome

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481