

XMU-MP-5

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	XMU-MP-5 acts as a selective inhibitor of ALK. It inhibits ALK-mutant Ba/F3 cells with an IC50 ranging from 4-50 nM and induces apoptosis in EML4-ALK Ba/F3 cells. Additionally, XMU-MP-5 demonstrates antitumor activity in mice.
Targets(IC50)	Apoptosis,ALK

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481