

TS-2

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	TS-2 is a fluorescent ligand for c-MycG4 with anticancer activity. It targets the G4 structure in the promoter of the c-Myc oncogene, inhibiting transcription of this cancer gene, which suppresses cancer cell growth and induces cell apoptosis (apoptosis).
Targets(IC50)	Apoptosis,Others,c-Myc,DNA/RNA Synthesis

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481