

LL-37 FK-13 TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	LL-37 FK-13 TFA, the TFA salt form of LL-37 FK-13, is an antimicrobial agent that inhibits <i>Trichomonas vaginalis</i> . It exhibits minimal hemolytic effects on human red blood cells and shows slight cytotoxicity towards human fibroblasts.
Targets(IC50)	Parasite

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481