

Estrogen receptor β /HDAC probe 1

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Estrogen receptor β /HDAC probe 1 (compound P1) is a near-infrared fluorescent probe designed for dual targeting of the Estrogen Receptor (Estrogen Receptor/ERR) β and the histone deacetylase HDAC.
Targets(IC50)	Estrogen Receptor/ERR,Others,HDAC

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481