

(S,R,R)-VH032-CO-C3-CONH-C2-PEG3-OH

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | (S,R,R)-VH032-CO-C3-CONH-C2-PEG3-OH is a conjugate of the E3 ligase ligand and linker for PROTACNCOA4 degrader-1 (E3LigaseLigand-Linker Conjugates), utilized in the synthesis of PROTACs. |
| Targets(IC50) | E3 Ligase Ligand-Linker Conjugates |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481