

SIYNFEKL TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | SIYNFEKL TFA, a variant of the primary MHC class I-restricted epitope SIINFEKL, serves as an antigenic peptide for stimulating specific T cells in experimental settings to study competitive interactions among T cells. It exhibits lower affinity for the OT-I T cell receptor (TCR), making it useful for detecting CD8+ T cells. |
| Targets(IC50) | Others,MHC |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481